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INDICATIONS

WINN'SAVE® devices for restraining or securing patients in a bed are subject to medical prescription. These products are designed and used for safety reasons for beds in the following cases:

- Physical restraint: to restrict the movements of an agitated and/or aggressive patient presenting cognitive or mental problems, to stop him/her falling and to protect staff or the patient from him/herself.
- Securing: protecting a confused patient from getting up unauthorised for medical reasons and for administering treatment.

PRECAUTIONS OF USE

This product must be attached by qualified personnel according to the manufacturer's instructions by following the safety regulations and procedures in line with the best practice recommendations applicable in health centres and regarding home care.

- Follow the instructions for use: restraining or securing in a chair depending on the indication, choice of size depending on waist circumference, adjustment of the length of the removable pelvic strap, comfortable tightness that does not restrict normal breathing at rest, specific surveillance of the skin in direct contact (rubbing, maceration).
- Before use, check that the characteristics of the restraining device are intact (eyelets, loops, seams, buckles, pelvic strap).
- Always attach the waist restraint to the pelvic strap using the side straps, unless contraindicated by the prescriber.
- In the event of direct contact with a patient in isolation for a skin condition (contact precautions), the device will be considered for use by said patient only, unless reuse is authorised after treatment by the hospital hygiene department.
- Use only the closing systems provided or approved by the manufacturer.

SPECIFIC SURVEILLANCE

- Prohibit any alteration or repair to the restraining or securing devices.
- Check that an agitated patient cannot access or handle the closing systems.
- Check the restraining devices used and stocked for the intention of physical restraint at regular intervals. They must be clean, complete, in good condition and kept such that they can be used swiftly.
- Keep unused devices out of patients' reach.
- Check that cleaning and maintenance are carried out in accordance with the instructions given in the user's manual.
- Check that the buckle for adjusting the length of the removable pelvic strap is protected between the 2 rear fabrics of the waist restraint.

Designation	Measurements	Closing system		
		Magnetic	Mechanic	Buckle and Velcro
Waist and pelvic restraint for the bed with side straps (S)	45-85 cm*	VCBTLT002-10	VCBTLT102-10	
Waist and pelvic restraint for the bed with side straps (M)	70-118 cm*	VCBTLT002-20	VCBTLT102-20	
Waist and pelvic restraint for the bed without side straps (M)	70-118 cm*	VCBTLT002-21	VCBTLT102-21	
Waist and pelvic restraint for the bed with side straps (L)	94-158 cm*	VCBTLT002-30	VCBTLT102-30	
Waist and pelvic restraint for the bed without side straps (L)	94-158 cm*	VCBTLT002-31	VCBTLT102-31	
Optional: bed rail for ankle trap	-	VCBTLT003-01	VCBTLT103-01	
Chest harness extension	-	VCBTLT004-01	VCBTLT104-01	
Chest harness restraint	20-50 cm	VCBTLT004-20	VCBTLT104-20	
Universal strap restraints for wrists and ankles (x2)	5-40 cm**	VCBTLT006-00	VCBTLT106-00	
One-size strap restraints for wrists (x2)	10-30 cm**	VCBTLT011-00	VCBTLT111-00	
Universal strap with buckle (x2)	5-40 cm**			VCNLT007-00
Waste strap (Buckles, Velcro) - (M)	65-85 cm*			VCNLT010-20
Waste strap (Buckles, Velcro) - (L)	85-120 cm*			VCNLT010-30

*Measurement corresponding to waist circumference

**Measurement corresponding to ankle and wrist circumference

Reference	Designation
252820004	Complete magnetic lock (key included)
252820006	Magnetic lock (without key)
302880012	Magnetic key
252820002	Mechanical lock (disc + pivot)
4099020014	Pen-key

INSTRUCTIONS

Closing

Put the metallic tip through the eyelet hole in the fabric strap. Insert the strap eyelet intended to close the restraining device. Close by clicking the button on the tip.


Precautions: when installing a restraining device, check that the closing system is correctly locked. To facilitate opening, release the pressure applied to the button.

Device for restraining or securing the patient's chest to the bed


Precautions: the waist restraint straps are fastened in the middle of the bed's long sides. Tighten until crushing the edge of the mattress. Lock the moving functions of the bed if possible (back and leg rests). Adjust the length of the removable pelvic strap before use according to the patient's morphology.

Chest restraint with or without extensions

Precautions: In the case of a medical bed with adjustable height, fasten the device to the mobile parts that move with the patient. In the case of a non-medical bed, fasten the device to the bed base or frame. Lock the moving functions of the bed if possible (back and leg rests).


Device for restraining or securing wrists

Secure wrists on the waist restraint strap or to the fixed side part of the bed.

Secure wrists using the buckle and fasten to the fixed structure of the bed.


Ankle restraints

Restrain ankles by fastening to the fixed structure of the bed rail. In the case of severe agitation, the bed rail can be secured using an extension secured to the middle of the rail and bed frame. Lock the moving functions of the bed base.

SAFETY LOCKS

MECHANICAL SAFETY LOCK

The mechanical safety lock is designed to prevent the patient can liberate himself, but can be opened simply by the medical personnel with a pen with a conical tip similar to the supplied with the closure.

The difference in diameters and shape of the holes of the disc prevents the pivot can be placed in the wrong hole.

CLOSE THE LOCK

1. Taking the closing disc in the hand, insert the center hole in the pivot, keeping the marked side at the top.
2. Press down until you hear "click" it indicates that the spring has been housed in the pivot slot.
3. Check, stretching upwards, that the disc has been properly secured.

OPEN THE LOCK

4. Insert the tip of the pen into the hole pointed by the arrow on one side of the disc.
5. Press against the resistance of the internal spring and pull the disc out.

IMPORTANT:

Always use a pencil similar to the pen supplied with the locks.

MAGNETIC SAFETY LOCK

The magnetic safety lock is designed to prevent the patient can free himself, but lets be opened simply by the medical personnel with the magnetic key supplied.

CLOSE THE LOCK


1. Taking the closing button in the hand, insert the center hole in the pivot.
2. Check, stretching upwards, that the button has been properly secured.

OPEN THE LOCK

3. Place the key provided on the top of the button.
4. Slide up proving that the magnet has liberated the button.

IMPORTANT: In patients with pacemakers should consider instructions' pacemaker manufacturer, regarding to possible damage because of strong magnets.

EMERGENCY OPENING: The button has two holes on the top. It allows an emergency opening by introducing, at the same time, a sharp object of appropriate diameter. (For example: a paperclip)

IMPORTANT SAFETY NOTES:

In case of blocking, cut the textile part of the restraining device, never cut the lock.

Warning: if locks are immersed in liquid, the proper functioning could be affected.

If incontinence, locks must be cleaned and disinfected immediately.

When washing, the anti-corrosive treatment of the locks can be seriously damaged. Therefore it is not recommended to wash the locks in a washing machine and certainly not at a temperature above 80 °C

We recommend cleaning the locks with any non-corrosive disinfectant product, like for example: a surface disinfectant with alcohol base <30%. Some trademarks are V5 Anioxo-spray or Bacterisan B1.

Closures must be removed before washing the restraint.

Never use the safety locks until you have read and understood this instructions.

The improper use of restraining devices could increase the risk of serious injury.

The restraining devices must be used ALWAYS under medical prescription.

ALL THE CLOSURES MUST BE REPLACED AFTER 3 YEARS (manufacturing date is printed in the lock)

Care instructions for the fabric products only

	<p>Machine washing up to 90°. Cooldown principle recommended</p>
	<p>Disinfectant in bleach (chlorine) compatible up to 5000 ppm or 0.5% active chlorine depending on best practice recommendations.</p>
	<p>Do not iron</p>
	<p>Not suitable for dry cleaning</p>
	<p>Tumble-dryer authorised at medium temperature</p>
<p>Fire Retardant Treatment </p>	<p>All of the fabric materials have undergone Fire Retardant Treatment</p>
	<p>Clean with Detergent-Disinfectant or Alcohol-based surface disinfectant (< 30%), with CE marking</p>

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