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INDICATIONS

WINN'SAVE® devices for restraining or securing patients in a chair are subject to medical prescription. These products are designed and used for safety reasons for armchairs and/or wheelchairs in the following cases:

- Physical restraint: to restrict the movements of an agitated and/or aggressive patient presenting cognitive or mental problems, to stop him/her falling and to enable treatment to be administered.
- Keeping patients upright: to guarantee the physical stability of a patient who is disorientated or has weak muscle tone and to increase his/her safety and comfort or administer treatment.

PRECAUTIONS OF USE

This product must be installed by qualified personnel according to the manufacturer's instructions by following the safety regulations and procedures in force in the health centre or regarding home care.

- Follow the instructions for use: restraining or securing in a chair depending on the indication, choice of size depending on waist circumference, adjustment of the length of the removable pelvic strap, comfortable tightness that does not restrict normal breathing at rest, specific surveillance of the skin in direct contact (rubbing, maceration).
- Before use, check that the characteristics of the restraining device are intact (eyelets, loops, seams, buckles, pelvic strap).
- Always attach the waist restraint to the pelvic strap.
- This product must not come into contact with injured skin.
- Use the closing systems provided or approved by the manufacturer.
- In the event of direct contact with a patient in isolation for a skin condition (contact precautions), the device will be considered for use by said patient only, unless reuse is authorized after treatment by the hospital hygiene department.

SURVEILLANCES SPECIFIQUES

SPECIFIC SURVEILLANCE

- Prohibit any alteration or repair to the restraining or securing devices.
- Place and secure the waist belt such that the patient is sitting right back in his/her seat, with the lumbar region against the chair back.
- Check that an agitated patient cannot access or handle the closing systems.
- Tilt the seat or chair back by 10° for patients with weak muscle tone.
- Adjust the chair foot rest correctly (height, angle, depth) so that the patient cannot slide forwards.
- If necessary, use extra devices like cushions, whose material limits sliding of upper body supports and hip abduction braces between the thighs or customized seats.
- Assess the merits of using restraining or securing devices on the arm/wheelchair daily.
- Check the restraining devices used and stocked for the intention of physical restraint at regular intervals. They must be clean, complete, in good condition and kept such that they can be used swiftly.
- Keep unused devices out of patients' reach.
- Check that cleaning and maintenance are carried out in accordance with the instructions given in the user's manual.
- Check that the buckle for adjusting the length of the removable pelvic strap is protected between the 2 rear fabrics of the waist restraint.

Designation	Measurements	Closing system		
		Magnetic	Mechanic	Sliding buckle
Optional: chest restraining harness for arm/wheelchair	One-size (U)	VCBTFT005-00		
Chest and pelvic restraint (S) for arm/wheelchair	45-85 cm*	VCBTFT005-10	VCBTFT105-10	-
Chest and pelvic restraint (M) for arm/wheelchair	70-118 cm*	VCBTFT005-20	VCBTFT105-20	-
Chest and pelvic restraint (L) for arm/wheelchair	94-158 cm*	VCBTFT005-30	VCBTFT105-30	-
Accessory - removable pelvic strap replacement	-	VCBTLF013-01		
Jacket with pelvic strap - (S)	50-80 cm	-	-	VCNTFT008-10
Jacket with pelvic strap - (M)	70-100 cm	-	-	VCNTFT008-20
Jacket with pelvic strap - (L)	90-130 cm	-	-	VCNTFT008-30
T-shaped waist restraint - (S)	50-80 cm	-	-	VNCTFT012-10
T-shaped waist restraint - (M)	70-100 cm	-	-	VNCTFT012-20
T-shaped waist restraint - (L)	90-130 cm	-	-	VNCTFT012-30

*Measurement corresponding to waist circumference (min./max.)

Reference	Designation	Reference	Designation
252820004	Complete magnetic lock (key included)	252820002	Mechanical lock (disc + pivot)
252820006	Magnetic lock (without key)	4099020014	Pen-key
302880012	Magnetic key		

INSTRUCTIONS

Magnetic closing


Put the metallic tip through the eyelet hole in the fabric strap. Insert the strap eyelet intended to close the restraining device. Close by clicking the button on the tip. The magnetic key will open the device by removing the button.

Precautions: when installing a restraining device, check that the closing system is correctly locked.
To facilitate opening, release the pressure applied to the button.


Device for restraining or securing the patient's waist and/or pelvis to the wheelchair

Waist strap + pelvic strap – Magnetic closing – Ref.: VCBTFT005-10 / 20 / 30

Place the waist strap with its perineal part preadjusted by making sure that you are not obstructing the patient's breathing when you tighten the waist strap.

With chest harness for restraint – Magnetic closing – Ref.: VCBTFT005-00

Place the chest strap with the perineal part then secure the front part of the harness first and then the rear part, making sure the fastening point is in the middle.

Jacket + pelvic strap – Sliding buckle and closing – Ref.: VCNTFT008-10 / 20 / 30


Fig. 1

Fig. 2

Fig. 3

Fig. 4

Install the jacket on the seat. (Fig. 1)

Sit the patient down. Lift the perineal part strap and each of the straps passing over the patient's shoulders to the back of the arm/wheelchair. (Fig.2)

Adjust the length using the sliding buckle. (Fig.3)

Secure everything using the quick fasteners. (Fig.4)

T-shaped waist and pelvic restraint – Sliding buckle and closing – Ref.: VCNTFT012-10 / 20 / 30

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Install the device on the arm/wheelchair. (Fig. 1)

Place the waist strap around the patient's waist and cross the straps behind him/her level with his/her lumbar region. (Fig. 2, 3)

Attach the straps on each side by adjusting the tightness – without compressing the patient's waist. (Fig. 4)

Fasten the end of the perineal part to the straps behind the arm/wheelchair. (Fig.5)

SAFETY LOCKS

MECHANICAL SAFETY LOCK

The mechanical safety lock is designed to prevent the patient can liberate himself, but can be opened simply by the medical personnel with a pen with a conical tip similar to the supplied with the closure.
The difference in diameters and shape of the holes of the disc prevents the pivot can be placed in the wrong hole.

CLOSE THE LOCK

1. Taking the closing disc in the hand, insert the center hole in the pivot, keeping the marked side at the top.
2. Press down until you hear "click" it indicates that the spring has been housed in the pivot slot.
3. Check, stretching upwards, that the disc has been properly secured.


OPEN THE LOCK

4. Insert the tip of the pen into the hole pointed by the arrow on one side of the disc.
5. Press against the resistance of the internal spring and pull the disc out.

IMPORTANT:

Always use a pencil similar to the pen supplied with the locks.

MAGNETIC SAFETY LOCK

The magnetic safety lock is designed to prevent the patient can free himself, but lets be opened simply by the medical personnel with the magnetic key supplied.

CLOSE THE LOCK


1. Taking the closing button in the hand, insert the center hole in the pivot.
2. Check, stretching upwards, that the button has been properly secured.

OPEN THE LOCK

3. Place the key provided on the top of the button.
4. Slide up proving that the magnet has liberated the button.

IMPORTANT: In patients with pacemakers should consider instructions' pacemaker manufacturer, regarding to possible damage because of strong magnets.

EMERGENCY OPENING: The button has two holes on the top. It allows an emergency opening by introducing, at the same time, a sharp object of appropriate diameter. (For example: a paperclip)

IMPORTANT SAFETY NOTES:

In case of blocking, cut the textile part of the restraining device, never cut the lock.

Warning: if locks are immersed in liquid, the proper functioning could be affected.

If incontinence, locks must be cleaned and disinfected immediately.

When washing, the anti-corrosive treatment of the locks can be seriously damaged. Therefore it is not recommended to wash the locks in a washing machine and certainly not at a temperature above 80 °C

We recommend cleaning the locks with any non-corrosive disinfectant product, like for example: a surface disinfectant with alcohol base <30%. Some trademarks are V5 Anioxy-spray or Bacterisan B1.

Closures must be removed before washing the restraint.

Never use the safety locks until you have read and understood this instructions.

The improper use of restraining devices could increase the risk of serious injury.

The restraining devices must be used ALWAYS under medical prescription.

ALL THE CLOSURES MUST BE REPLACED AFTER 3 YEARS (manufacturing date is printed in the lock)

Care instructions for the fabric products only

	<p>Machine washing up to 90°. Cooldown principle recommended</p>
	<p>Disinfectant in bleach (chlorine) compatible up to 5000 ppm or 0.5% active chlorine depending on best practice recommendations.</p>
	<p>Do not iron</p>
	<p>Not suitable for dry cleaning</p>
	<p>Tumble-dryer authorised at medium temperature</p>
<p>Fire Retardant Treatment </p>	<p>All of the fabric materials have undergone Fire Retardant Treatment</p>
	<p>Clean with Detergent-Disinfectant or Alcohol-based surface disinfectant (< 30%), with CE marking</p>

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