
WINCARE

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
DEVICE INTENDED USE

This high resilience foam static mattress is intended to be used in combination with a medical bed for specific pressure ulcer prevention purposes due to a loss of autonomy and/or deterioration of health resulting in reduced mobility and the risk of pressure ulcers.

INDICATIONS

Bedridden patients with reduced mobility and "Low to Medium" risk of pressure ulcers (according to the Braden scale or any other validated scale and on medical advice). Assists in the treatment of stage 1 to 2 pressure sores according to EPUAP, up to stage 4 outside the support area.

CONTRAINDICATIONS

Patient weight level below or above the limits defined in the technical data table. For use on medical imaging table.

ADVERSE SIDE EFFECTS

Not documented.

Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the device poses a serious risk or has been tampered with.

MEDICAL DEVICE COMPOSITION/TECHNICAL DATA

| References Models and versions* | VA103MB14 | VA103MB14RC | VA103MB14/HIP | VA103MB14/HIPH |
|---------------------------------|---------------------------------------|-------------|---------------------------------------|----------------|
| Sizes | 195x85x14 cm | | 195x85x14 cm | |
| Lower specs | HR: 42.9 kg/m ³ 4.4 kPa | | HR: 42.9 kg/m ³ 4.4 kPa | |
| Patient weight min and max | 40 to 120 kg | | 40 to 120 kg | |
| Foam material | Dermalon FR (Fire resistant) | | 100% coated polyurethane on PES weave | |

The references mentioned can be completed with a suffix to constitute a commercial reference (example VA103MB14/TH)

CLINICAL BENEFIT, DEVICE PERFORMANCE, MECHANISM OF ACTION

DEVICE PERFORMANCE CHARACTERISTICS: By immersion, this support relieves pressure on the body by improving blood circulation and thus reducing the risk of pressure ulcers on the supporting areas of the body.

EXPECTED CLINICAL BENEFITS: Maintenance of tissue oxygenation at the cutaneous and subcutaneous tissue level in contact with the support.

INFORMATION FOR HEALTHCARE PROFESSIONALS: Monitor the condition of the patient's skin in contact with the support several times a day. Use pressure relief devices or positioning systems in patients with pressure ulcer(s).

PREREQUISITES BEFORE USE AND INSTRUCTIONS FOR USE

TRAINING AND QUALIFICATION OF THE DEVICE USER: user training must be carried out by individuals trained and approved by the economic operators concerned, notably in terms of safety and non-compliance reporting.

DEVICE installation: The mattress is delivered with its draw sheet or removable zipped cover. The product is ready to be installed. The surface of the foam composed of studs ("alveolar") must be positioned in contact with the patient's body. The protection must be dry.

Preventive maintenance: Carry out a regular visual inspection of the foam condition: visible sagging of material and a non-consistent and very slow return of the foam shape are ageing criteria that compromise the properties of the support. Annual inspection of the condition of the cover (surface aspect and zippers) by exposing the inner surface to a light source: no holes and/or tears. Replace the cover if its surface appearance changes.

CLEANING AND DISINFECTION

| Promust cover PU and HDPV |  |  |  |  |  |  |
|------------------------------|---|---|---|---|---|--|
| Moderate washing up to 90°C | Maximum chlorine concentration of 5000ppm | Do not iron | Do not dry clean | Dry with reduced thermal restriction | Use an approved surface disinfectant/detergent | |
| Dermalon draw sheet |  |  |  |  |  |  |

INFORMATION ON APPROPRIATE PROCESSES TO ALLOW ITS REUSE

The product must exhibit physical and bacteriological cleanliness. "Promust" protection treated with a biocide substance at no risk to users. Avoid wet contact with the foam. Prohibit all scouring, stripping products or solvents and "piercing-sharp" objects in direct contact with the cover.

WARNINGS, PRECAUTIONS FOR USE, MEASURES REQUIRED**PRECAUTIONS FOR USE**

A bedsore is a lesion of varying depth in the skin, associated with excessive and prolonged compression of tissue between the body and a support surface. This excessive pressure can prevent blood circulation and lead to pressure ulcers, a dangerous lesion because it develops from the inside of the body to the outside, and is therefore initially invisible. A pressure sore can have several forms: a simple persistent rash for more than a day, hardening of the skin, a wound of varying depth that can in severe cases affect muscles or underlying bone.

Pressure (max) in mmHg before ageing, at the maximum validated weight, measured flat (0°)

| SUITABILITY FOR USE | Max weight User | Head | Buttocks | Heel |
|---------------------------|-----------------|------|----------|------|
| APLOT mattress (mm of hg) | 140 kg | 54 | 46 | 92 |

WARNING

(Re)assess the risk of patients being trapped in the non-moving parts of the medical bed associated with the "therapeutic mattress and accessories and the articulated bed base positions". The use of segmented or full length side barriers must be required. External cardiac massage not compatible without a board placed between the thorax and the top surface of the mattress

STORAGE, HANDLING, DISPOSAL**STORAGE AND USAGE CONDITIONS**

The supports must ideally be stored flat, away from direct light and excessive humidity.

| | Use | Storage |
|----------------------|---|---------------------------------|
| Temperature range | +15°C to +45°C / +59°F to +113°F | -25°C to 70°C / +13°F to +158°F |
| Humidity range | 30% - 70% | 30% - 95% |
| Atmospheric pressure | 50kPa - 106 kPa // 7.25 Psi / 15.37 Psi | |

ROLLED MATTRESSES

APLOT: 30 mn → 2h : This time may be increased if the storage temperature is lower than 15°C and/or in line with the storage period (0 to 5 months). The return to a normal visual appearance of the mattress under normal usage temperature conditions may take up to 24 hours for ALOVA, APLOT and MORFEA mattresses but this does not present any risk for their use. Small creases may be visible on the lower surface of the CIC cover after unrolling. They present no risk to the patient and do not affect the efficiency of the final product.

SERVICE LIFE

Highly resilient mattresses have an estimated service life of 6 years.

The rolled-compressed mattress life span is 6 months (see date on package label).

PRODUCT DISPOSAL

Do not dispose of the product outside of the dedicated area. Respect the established recycling systems in your country.

GUARANTEE

APLOT mattresses have a 2-year guarantee. This guarantee starts from the date of purchase of the product from your distributor. This guarantee does not cover normal wear and tear of the product and its protection. Please contact your distributor about a defective product. They will carry out the necessary procedures with our company to proceed with either repair or exchange.


Class I medical device according to Regulation (EU) 2017/745.

Management systems: Quality: ISO 13485:2016 - Environmental: ISO 14001:2015

Technical compliance certificates (FCBA): ESC N° 16-005 of 10 June 2016

Contents of package

1 mattress
1 protective cover
1 user manual

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