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WIN CARE

User manual Viscoelastic foam cushions

ALOVA CURVED – ALOVA ERGONOMIC **DEVICE INTENDED USE**

Pressure sore prevention in prolonged sitting position.

INDICATIONS

A seated patient with reduced mobility, no support asymmetry and no risk of slipping forward, with a "Moderate to High" pressure ulcer risk (according to the Braden or other approved scale and in accordance with medical opinion). Trunk and/or lower limb paralysis, history of ischial or sacrococcygeal pressure ulcer, lower limb tenderness, spasticity, vascular amputation.

CONTRAINDICATIONS

Patient with sagittal or frontal stability disorders. Sacrococcygeal and/or ischial pressure ulcers. Patient weight levels above the limits set in the technical data table.

ADVERSE SIDE EFFECTS

Trunk dynamic instability. Pelvic slippage in the sagittal plane. Loss of efficiency in "sitting" transfers. Lateral trunk transfer instability.

Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the device poses a serious risk or has been tampered with.

MEDICAL DEVICE COMPOSITION/TECHNICAL DATA

| | VCAG | | | | VCAE |
|----------------------------|---|--|-------------------------------------|--|--------------|
| Versions | VCAG3636, CAG3841, VCAG4041, VCAG4340, VCAG4341, VCAG4043, VCAG4343, VCAG4347, VCAG4540, VCAG4544, VCAG5044 | VCAG5250, VCAG5450, VCAG5650, VCAG6050 | VCAG4343XL, VCAG5044XL | VCAG3841(BP), VCAG4041(BP), VCAG4544(BP), VCAG4347(BP), VCAG5044BP, VCAG4341(BP), VCAG4343(BP) | VCAE/M4343 |
| Dimensions | The Length and width dimensions are contained in the reference | | | | |
| Thickness | 6 cm | 7 cm | 9 cm | 6 cm | 9 cm |
| Surface | Viscoelastic > 75 kg/m ³ | | | | |
| Foam base | - | - | HR > 34 kg/m ³ | - | - |
| Protection | Removable cover in PU 30% / Polyester70% and non-slip lower side | | | | Promust PU |
| Patient weight min and max | 30 to 120 kg according to reference | 120 to 200 kg according to reference | 30 to 135 kg according to reference | 30 to 135 kg according to reference | 30 to 100 kg |

The references mentioned can be completed with a suffix to constitute a commercial reference (example VCAG4343-TH)

CLINICAL BENEFIT, DEVICE PERFORMANCE, MECHANISM OF ACTION

DEVICE PERFORMANCE CHARACTERISTICS: The viscoelastic foam material allows the pelvis to be immersed and the areas opposite bone protuberances (ischium, coccyx, trochanter) to be enveloped.

EXPECTED CLINICAL BENEFITS: Maintenance of tissue oxygenation in the cutaneous and subcutaneous tissues in contact with the support.

INFORMATION FOR HEALTHCARE PROFESSIONALS: Monitor the condition of the patient's skin in contact with the cushion several times a day. Do not seat patients with ischial and coccygeal pressure ulcers. Check that the seat of the chair is in good condition. In the case of a wheelchair, check the height adjustment of the armrests, footrests, (...) and performance in terms of propulsion.

PREREQUISITES BEFORE USE AND INSTRUCTIONS FOR USE

TRAINING AND QUALIFICATION OF THE DEVICE USER: user training must be carried out by individuals trained and approved by the economic operators concerned, notably in terms of safety and non-compliance reporting.

DEVICE installation: The cushion is delivered with its cover fitted. The product is ready to be installed. The upper side of the cushion in contact with the patient is made up of viscoelastic foam, the ischial curves and coccygeal discharge towards the rear and external rotator towards the front.

Do not use the ALOVA CURVE and ERGONOMIC cushion without the protective cover.

Preventive maintenance: Carry out a regular visual check of the condition of the cover (surface appearance and zippers) by exposing the inner surface to a light source: no holes and/or tears) and annually check the surface condition of the foam. Replace the protective cover if its surface appearance changes. The non-slip surface is in contact with the seat of the chair. The handle (located underneath the product) must face the front of the cushion.

CLEANING AND DISINFECTION

| Procedure |  |  |  |  |  |  |
|------------------|---|---|---|---|---|--|
| Cover Promust PU | Moderate washing up to 90°C | Maximum chlorine concentration of 5000ppm | Do not iron | Do not dry clean | Dry with reduced thermal restriction | Use an approved surface disinfectant/detergent |

INFORMATION ON APPROPRIATE PROCESSES TO ALLOW ITS REUSE

Before re-using for another patient, the product must exhibit physical and bacteriological cleanliness after application of the above procedures. Item treated with a biocide substance at no risk to users. www.winnocare.fr

Avoid direct damp cushion contact with the skin. Always protect the cushions with the supplied protective cover. Prohibit all scouring, stripping products or solvents and "piercing-sharp" objects in direct contact with the cover.

WARNINGS, PRECAUTIONS FOR USE, MEASURES REQUIRED**PRECAUTIONS FOR USE**

Patient with bone and/or muscle and/or skin trauma to the pelvic region and thigh in contact. Seating instability.

A bedsore is a lesion of varying depth on the skin, associated with excessive and prolonged compression of tissue between the body and a support surface. This excessive pressure can prevent blood circulation and lead to pressure ulcers, a dangerous lesion because it develops from the inside of the body to the outside, and is therefore initially invisible. A pressure sore can have several forms: a simple persistent rash for more than a day, hardening of the skin, a wound of varying depth that can in severe cases affect muscles or underlying bone.

Pressure (max) in mmHg before ageing, at the maximum validated weight, measured flat (0°)

| SUITABILITY FOR USE | User weight | Head | Ischions | Heel |
|--------------------------------|-------------|------|----------|------|
| ALOVA CURVE cushion (mm of Hg) | 100 kg | NA | 91.20 | NA |
| ERGONOMIC cushion (mm of Hg) | 100 kg | NA | 85.32 | NA |

STORAGE, HANDLING, DISPOSAL**STORAGE AND USAGE CONDITIONS**

The supports must ideally be stored flat, away from direct light and excessive humidity.

| | Use | Storage |
|-------------------|-----------------------------------|--------------------------------|
| Temperature range | +15°C to +45°C +59°F to +113°F | -5°C – 60°C +41°F to +140°F |
| Humidity range | 30% - 70% | 10% - 90% |

SERVICE LIFE

The cushion has an estimated service life of 6 years.

PRODUCT DISPOSAL

Do not dispose of the product outside of the dedicated area. Respect the established recycling systems in your country.

GUARANTEE

The guarantee duration for cushions is 3 years. This guarantee starts from the date of purchase of the product from your distributor. In the case of single patient use, the guarantee relates to the support and its protective cover.

This guarantee does not cover normal wear and tear of the product and its protective cover and is not a substitute for legal guarantees. Please contact your distributor about a defective product. They will carry out the necessary procedures with our company to proceed with either repair or exchange.

Technical compliance certificates

| FCBA - ESC NO. 19-006 | FCBA - ESC NO. 19-008 | FCBA - ESC NO. 13-007 | FCBA - ESC NO. 15-003 |
|---|--------------------------------------|--------------------------------------|----------------------------------|
| ALOVA Curved Standard Cushion ALOVA Curved Cushion Standard < 38 cm ALOVA Cushion Curved Standard > 45 cm | ALOVA Cushion Curved Technical XL | ALOVA Cushion Curved BP Technique | ALOVA Ergonomic Standard Cushion |

Class I medical device according to Regulation (EU) 2017/745.

Management system: Quality: ISO 13485:2016 - Environmental: ISO 14001:2015

Contents of package

- 1 cushion
- 2 covers. For the BP reference a cover with a housing for the pelvic end-piece.
- 1 user manual

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| Processo | | | | | | |
|------------------|-----------------------|---|-------------|------------------|---|---|
| Cover Promust PU | Moderate wash at 90°C | Maximum authorised chlorine concentration 5000ppm | Non stirare | Do not dry clean | Drying with reduced thermal constraints | Use of an authorised surface detergent/disinfectant product |

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(max) Pressure in mmHg before ageing, using the maximum validated weight, measured flat (0°)

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This warranty does not cover normal wear and tear and is not a substitution of legal guarantees. Please contact your distributor and present the incriminating product. The distributor will carry out the necessary procedure with our company in order to either repair or replace the product.

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