

4 Le Pas du Château
85670 SAINT PAUL MONT PENIT
TEL : +33 (0)2 51 98 55 64
FAX : +33 (0)2 51 98 59 07
Email : info@medicatlantic.fr
Site Internet : <http://www.winnocare.fr>

INSTRUCTION MANUAL

EVILENCE BED TABLE T3800



SUMMARY

- ① TRANSPORT
- ② STORAGE
- ③ ASSEMBLY
- ④ USE
- ⑤ MAINTENANCE
- ⑥ DISPOSAL

YOU MUST READ THIS BEFORE USING THE PRODUCT.

COPY TO BE KEPT.



WARNING

These instructions are part and parcel of the product, and they contain technical requirements which must not be reproduced fully or partially, or divulged or used without authorisation for competitive purposes, or given to third parties.

*Although all the information contained in these instructions have been carefully checked and can be taken as reliable, **MEDICATLANTIC** assumes no responsibility for any errors, omissions or misprints.*

① TRANSPORT

When transporting the table, it must be in a low position and packaged in cardboard. **IMPORTANT:** it is strictly forbidden to stack packages on the table.

② STORAGE

The table must be stored away from the damp in a dry room at 5–40°C.

③ ASSEMBLY

Non applicable


④ USE

The table is designed for eating and reading at.

THE LOAD ON ITS SURFACE SHOULD NEVER EXCEED 15kg.

4.1 GENERAL DESCRIPTION

The table is made up of a base (1), a leg (23), a bracket (2), a three-side rack (3), a large desktop surface (5), a small desktop surface (6) and a book-rest (7).



4.2 USE

To raise or lower the table, lift the lever (29) and adjust the top to the desired level.

4.3 PRECAUTIONS

The table should be used in accordance with the advice given in this document, and also making sure that no object gets in the way of its zone of movement.

Immediately wipe the surface if water gets on it, for example through a glass of liquid being spilt. Such water damage is not covered by warranty.

④ MAINTENANCE

5.1 TECHNICAL MAINTENANCE

- The table should be carefully examined at least once a year to check parts such as the wheels, or the height adjustment, and also to make sure the desktop and rail are firmly attached.

5.2 CLEANING

- Clean the structure with soapy water and a soft cloth or sponge, then rinse and wipe.
- All harsh products such as detergents, powders, solvents and bleach are banned.

5.3 DISINFECTION

- Use a disinfectant which kills bacteria, fungi and viruses, either by spraying it evenly over the surfaces or by applying it with a single-use cloth, or by spraying a disinfectant aerosol from 30cm (1ft) away.
- Leave to dry, and isolate the disinfected equipment from other non-disinfected equipment using a film indicating the date of disinfection. (IF RENTED)


IMPORTANT:

Adhere to the usage precautions given on the disinfectant products.

5.4 WARRANTY

- All our products are guaranteed not to have any manufacturing defects under normal conditions of use and maintenance.
- Not covered are man-hours corresponding to changing structures or parts under warranty.
- Regarding the warranty duration specific to each product, please check the general conditions of sale.
- In any correspondence regarding possible maintenance, it is essential to quote the identification information on the shower-trolley and its jack, if applicable.
- Replacement is carried out via the supplying of original parts within the duration of the warranty by our resale network determining the start of the period of warranty.
- In order to allow proper application of this warranty and also to avoid any billing, it is compulsory to return any faulty parts.

5.5 IDENTIFICATION



© DISPOSAL

The product must be disposed of if essential requirements are not met, i.e. when the product no longer has its initial characteristics and has not been recalled during manufacture. The types of products used to make the do not directly affect the environment; however, it is recommended you take steps to make the chair no longer useable for its intended purpose.