


# GELTOP® / GELSCAR®



Geltop cushion




Gelscar cushion

## GELSCAR CUSHION

### KEY FEATURES

- ▶ Viscofluid gel cushion in an isothermal neoprene cover
- ▶ 4 cords sewn on the 4 corners
- ▶ Patient weight: 50 to 110 Kg (7.8 to 17.3 st)
- ▶ 2 sizes available:  
40 x 40 x 3.5 cm (15.7 x 15.7 x 1.1 inches)  
and 43 x 45 x 3.5 cm (16.9 x 17.7 x 1.1 inches)
- ▶ Hole version available:  
dim. 40 x 40 x 3.5 cm (15.7 x 15.7 x 1.1 inches)



## GELTOP CUSHION

### KEY FEATURES

- ▶ Viscofluid gel cushion combined with a High Resilient foam base
- ▶ Semi-permeable cover fitted with a non-slip base
- ▶ Patient weight: 40 to 95 Kg (6.2 to 14.9 st)
- ▶ 2 sizes available: 40 x 40 x 5.5 cm (15.7 x 15.7 x 1.9 inches)  
43 x 45 x 5.5 cm (16.9 x 17.7 x 1.9 inches)
- ▶ Thickness: 1 cm of viscofluid gel and 4.5 cm of High-Resilient foam.



For homecare use: delivered with 2 covers  
For nursing homes use: delivered with 1 cover

## INDICATIONS

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT LOW TO MEDIUM RISK OF PRESSURE ULCER DEVELOPMENT

- ▶ Cushion made up of gel (GELSCAR) and gel/foam (GELTOP)
- ▶ Viscofluid gel properties allow reducing shearing forces applied on the skin
- ▶ 4 cords (GELSCAR) or non-slip base (GELTOP) to stabilize the cushion on the (wheel)chair

DESCRIPTION	REFERENCE	DIMENSIONS	PATIENT WEIGHT	FIRE STANDARDS
GELSCAR for a chair seat 38/42 wide and 39/42 deep.	VGLS01C	40 x 40 x 3.5 cm	50 to 110 kg	EN 591-1 and 2
PIERCED GELSCAR for a chair seat 38/42 wide and 39/42 deep.	VGLS03C	40 x 40 x 3.5 cm		
GELSCAR for a chair seat 40/45 wide and 44/47 deep.	VGLS01C/4345	43 x 45 x 3.5 cm	40 to 95 kg	
GELTOP for a chair seat 40/42 wide and 40/42 deep, cover in Techmapolar.	VAB40	40 x 40 x 5.5 cm		
GELTOP for a chair seat 43/45 wide and 43/47 deep, cover in Techmapolar.	VAB4345	43 x 45 x 5.5 cm		

2 years warranty

Cleaning instructions for covers:



Photos not contractually binding

FP000043/01\_2018



WINCARE Group - 200 rue Charles Tellier - Actiparc de Grézan - 30034 Nîmes Cedex 1 France  
Ph : +33 (0)4 66 02 15 15 - Fax : +33 (0)4 66 02 15 00 - Email : contact.askle@winnicare.fr - www.winnicare.fr  
In the interest of constantly improving its products, the WINCARE Group reserves the rights to modify the technical specifications of any product without prior notice  
Class I medical device in accordance with European Directive 93/42/EEC

