

# Winn'Motion 175<sup>®</sup>



## FOR PATIENTS WEIGHT UP TO 175 KG

### SAFETY

- ▶ Spreader bar padded and secured by suspension points fitted with a non-return system.
- ▶ Control box with emergency stop.
- ▶ Audible signal when battery is being charged.
- ▶ Electric emergency lowering system.
- ▶ Anti-collapse system.
- ▶ Jack can be released using an ergonomic rotary knob.
- ▶ Battery autonomy 40 cycles (24 volts) (1 cycle = 1 lift/1 lowering), time required for a full charge: 8 to 10 hours.

### CONVENIENT

- ▶ Light and compact for use at home or in collective care setting.
- ▶ Easy to dismantle into 2 parts.
- ▶ Easy to transport and store away.
- ▶ Easy to handle: twin-wheel castors Ø 75 mm at the front and Ø100 mm at the rear.
- ▶ Weight: 36 kg.
- ▶ 360° spreader bar rotation.
- ▶ Design castors (optional).

**Available slings:**  
Bath, Fletty sling, 3D comfort, with or without head support, 2 sizes (M or L)

#### FLETTY SLING WITH HEAD SUPPORT



#### BATH SLING WITHOUT HEAD SUPPORT




Emergency stop button



Secured 4-point spreader bar



Design castors



CONFORMS TO STANDARD EN ISO 10535: 2007 APPLICABLE TO HOISTS FOR TRANSFERRING PEOPLE WHO ARE DISABLED.

### SAFETY

- ▶ Emergency stop button on the side of the control box.
- ▶ Hoist is suitable for patients weighing up to 175 kg.
- ▶ Audible signal when the battery is running low.

### COMFORT

- ▶ Jack with anti-collapse system.
- ▶ Jack can be manually released using an ergonomic rotary knob.
- ▶ Feet can be manually moved further apart using a pedal system.

DESCRIPTION	A	B	C	D	E	F	MAX PATIENT WEIGHT	FR CASTORS	RE CASTORS	STANDARDS	MOTORIZATION	WEIGHT
Winn'Motion 175	59 cm	185 cm	110 cm	65 cm	84 cm	12 cm	175 kg	7.5 cm	10 cm	ISO 10535	LINAK IP 54	36 kg

Hoist colors: RAL 7037 grey  
Remote control, electric and mechanical parts: 5-year warranty  
Batteries, slings: 1 year

Photos not contractually binding

FP000062/01.2018



WINNCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France  
Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatlantique.fr - www.winnicare.fr  
In the interest of constantly improving its products, the WINNCARE Group reserves the rights to modify the technical specifications of any product without prior notice  
Class I medical device in accordance with European Directive 93/42/EEC

