

ALOVA™ XXL



Braden scale score



KEY FEATURES

- ▶ Single-piece mattress made up of:
 - 7 cm molded visco elastic top layer: density: $\geq 75 \text{ Kg/m}^3$; hardness: 1.7 kPa
 - 10 cm High Resilient (HR) foam base: density: 42 Kg/m^3 ; hardness: 4.5 kPa
- ▶ Patient weight: 135 to 270 Kg
- ▶ Dimensions: VMA120BAR: 195 x 117 x 17 cm (76.7 x 46 x 6.6 inches)
VMA140BAR: 198 x 137 x 17 cm (77.9 x 53.9 x 6.6 inches)
- ▶ Mattress weight: 29.5 Kg (VMA120BAR) and 38.9 Kg (VMA140BAR)
- ▶ Cover available: Promust PU HD

KEY BENEFITS

- ▶ Visco elastic foam provides excellent pressure re-distribution properties: homogeneous immersion of the body into the foam, thus increasing the bearing surface and reducing pressure applied in areas at risk of pressure ulcers.
- ▶ Suitable for large width sleeping surfaces of 120 or 140 cm width medical beds for overweight or obese persons.
- ▶ Visco elastic foam provides patient comfort .

Note: the medical prescription must include the weight of the person.

INDICATIONS

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT MEDIUM TO HIGH RISK OF PRESSURE ULCER DEVELOPMENT, BEDRIDDEN MORE THAN 15 HOURS A DAY.

HELP IN TREATING EXISTING PRESSURE ULCERS FROM STADE 1 TO 2 COMBINED WITH TECHNICAL AID DEVICES.

PATIENT WEIGHT: 135 TO 270 KG/ 21.2 TO 42.5 ST

- ▶ Single-piece mattress made up of visco elastic foam top layer on High Resilient (HR) polyurethane foam
- ▶ ALOVA XXL: therapy and comfort in large width bed
- ▶ Removable and washable integral cover

DESCRIPTION	ITEM REFERENCES ACCORDING TO PACKAGING CARDBOX	FIRE STANDARDS
ALOVA™ XXL 120 WIDTH, PROMUST PU HD COVER	VMA120BAR	EN 597-1 and 2 GPEM D1 90 classe D
ALOVA™ XXL 140 WIDTH, PROMUST PU HD COVER	VMA140BAR	

3 years warranty

Cleaning instructions: Promust PU HD cover

FP000037/08.2018

Photos not contractually binding

