



BED ACCESSORIES

1



2



3



4



5



6



7



8




9



10




11



12



13





"Medical device accessory means any item which, without itself being a medical device, is intended by its manufacturer to be used with a given medical device or devices to enable the use of the medical device or devices in accordance with its or their intended purpose or to contribute specifically and directly to the medical function of the medical device or devices in accordance with its or their intended purpose." according to Regulation (EU) 2017/745 on medical devices (art. 2).

DESCRIPTION		ALDRYS	X'PRESS (LOW)	X'PRIM	MÉDICALYS	AERYS	XXL DIVISYS	XXL X'PRESS	DUO DIVISYS	KALIN	DIMENSIONS	FINISH/ MATERIALS
1	Epoxy Side Rail A639-00	•	•	•	•	•	•	•	•	-	Length in unfolded position: 1.47 m Height in unfolded position: 37.5 cm in relation to the bed base	Epoxy
2	Footboard rail A671-00	•	•	•	•	•	•	•	•	-	Length: 70 cm Height: 37.5 cm from the base	Epoxy
3	Wooden Side Rail A676-00	•	•	•	-	•	•	•	•	Specific Model	Total length of the bed Height: 37.5 cm from the base	Polymer Coating Chêne Lindberg or lacquer
4	Aluminium Half Rail A665-00/A666-00	•	•	•	•	•	•	•	•	-	3 heights of use	Aluminium/ Solid wood, lacquered or varnished
5	Urinal Holder A260-00	•	•	•	•	•	•	•	•	-	-	Epoxy
6	Reading Lamp 61300158 / 61300159	•	•	•	•	•	•	•	•	•	Mounted on the rod holder	Epoxy
7	Telescopic IV Stand with 2 or 4 hooks A1700 / A170001	•	•	•	•	•	•	•	•	•	2 hooks: fixed height 120 cm 4 hooks: adjustable height 115-120 cm	Stainless Steel
8	IV Stand on Castors with 2 or 4 hooks A8400 / A840001	•	•	•	•	•	•	•	•	•	Adjustable height: 72 to 141cm	Plastic Hooks
9	Flexible Remote Control Holder A230-00	•	•	•	•	•	•	•	•	-	-	Magnet / Epoxy
10	Urine Bag Hook A5800	-	-	-	•	-	-	-	-	-	-	Stainless Steel
11	Wall Bumper	A675-00	-	A551-00	A633	A636-00 (AÉRYS 90/100) A663-00 (AERYS XXL) A670-00 (AÉRYS without lateral movement)	A551-00	-	-	Standard	-	Epoxy
	Adjustable Wall Bumper A684-00	•	-	•	-	-	•	-	-	•	-	Epoxy
12	Bed Base Extension	•	•	•	•	•	•	•	-	-	Depth: 20 cm	Epoxy
13	Aerys Base Extension	-	-	-	-	90 cm wide bed only	-	-	-	-	Depth: 20 cm Length: 90 cm	Epoxy
14	Support Handle A575-00	•	•	•	•	•	•	•	•	-	Height: 35 cm Length: 29 cm	Epoxy
15	Angled Lifting Pole A622-00	•	•	•	•	•	•	•	•	Specific Model	-	Epoxy

